



1.2 BIOTECTUM

NEW GENERATION OF LAMINAR FLOW CABINETS

Class II microbiological safety cabinet, Biohazard,
made according to the EN12469 norm



made in
POLAND



BIOTECTUM





ecoMove

Energy demand reduced by more than 78 %, especially quiet operation, below 48 dB



Ergonomics and safety

Ergonomically designed control panel, iMove knob, intelligent system of warning and self-diagnosis



Equipment

Standard and optional



Effectiveness

Air flow velocities, structure



Technical data

Type of motors, power consumption, dimensions, illumination



AutoProtect, sterilization

Emergency power system, automatic and semi-automatic sterilization

05

07

08

09

10

11

BIOTECTUM 1.2

New generation of laminar flow cabinets

light bar
in LED technology,
visual information
on the operating state

media
connections
on the top
of the cabinet



3 EC type
motors

intelligent
management
system
iMove

electrically
adjustable
front glass

valves
and sockets
inside
the cabinet

cabinet interior made of
stainless steel of class 0H18
(DIN 1.4301)
in seamless construction
with rounded sides


worktop and
the lower bowl
made entirely of
stainless steel
of class 0H18
(DIN 1.4301)

V-shaped
air inlet



CLASS II MICROBIOLOGICAL SAFETY CABINET, BIOHAZARD, ACCORDING TO THE EN12469 STANDARD, WITH VERTICAL LAMINAR FLOW

Energy demand reduced by more than 78 %



A power consumption [W] ¹

B power consumption on an annual basis [kWh] ²

C charges for electricity on an annual basis [PLN] ³

D charges for electricity in terms of 10 years [PLN] ³, reduction of electricity costs by an amount of ~1000 EUR

1) In the safe work mode according to EN 12469 norm


2) Based on operation for 8 h/day, 5 days during the week, 52 weeks.

An additional assumption has been made for BioTectum, that 20 % of cabinet operation will be realized in the SUSPEND mode

3) average cost of 1 kWh = PLN 0.55, data source:

http://ure.gov.pl/ftp/ure-kalkulator/ure/formularz_kalkulator_html.php

CO₂ emissions reduced by more than 75 %





ecoMove TECHNOLOGY

During production of the new generation laminar flow cabinets a special emphasis was placed on improving the work safety, comfort and maximum reduction of energy demand. Thanks to applying three electronically commutated motors and a precise control, we gained an environment- and user-friendly product (low power consumption, low noise level and heat emission).

- 3 EC-type motors (electronically commutated)
- low power consumption < 120 W
- low noise level < 48 dB
- increased energy savings SUSPEND mode
- auto compensation of air flow velocity
- intelligent, dedicated, microprocessor based control system
- rigid construction based on an aluminum frame





ALARM SYSTEMS

The maximum safety is provided by acoustic and visual alarm messages with the possibility of approving the acoustic alerts by the operator of the laminar cabinet.

TYPES OF ALARM MESSAGES

- exceeded working time of UV lamp
- exceeded working time of filters
- exceeded working time without service inspection
- glass outside the working position
- exceeded air temperature in the cabinet
- interrupted sterilization process
- BEZPIECZNA mode inactive
- detected power outage
- operation on the UPS
- fans working on the peak performance
- barrier air velocity is too low

- air velocity in the laminar curtain is too low
- exceeded maximum air velocity in the laminar curtain

SELF-DIAGNOSIS SYSTEMS

Fast verification of the failure enables rapid response service.

TYPES OF DIAGNOSTIC MESSAGES

- sensor failure - air velocity in the laminar curtain
- sensor failure - barrier air velocity
- temperature sensor failure
- glass failure w kurtynie laminarnej

The type of displayed messages depends on the equipment options of the laminar cabinet.

EQUIPMENT / VERSION

COMFORT ADVANTAGE AUTOPROTECT


ull worktop, non-sectioned, removable	●	●	●
stylistic elements in orange RAL 2000	●	●	●
permanently mounted UV lamp	●	●	●
driver equipped with decontamination program	●	●	●
3 electrical sockets (2 on the left, 1 on the right-hand side of working chamber)	●	●	●
port for the DOP test of HEPA filters	●	●	●
MODBUS RTU communication protocol	●	●	●
cooperation with BMS (Building Management System)	●	●	●
RS-485 port USB port	●	●	●
clock/date	●	●	●
adjustment of the display contrast	●	●	●
frame for the cabinet made of closed profiles	○	●	●
valve for combustible gases	○	●	●
valve for technical gases	○	●	●
temperature sensor in the working chamber	○	●	●
adjustable illumination of the working area	○	●	●
AutoProtect mode	○	○	●
worktop, removable, sectioned	○	○	○
valve for combustible gases with solenoid valve (activated from the control panel)	○	○	○
valve for technical gases with solenoid valve (activated from the control panel)	○	○	○
control of an external exhaust fan	○	○	○
stylistic elements in other colour	○	○	○

● standard ○ option




Linear speed of the air in the working area (under the supply filter)

the test was made 0,10 m above the glass edge in 12 points
in the horizontal plane *



Linear speed of the air barrier (access window)

Velocity in reduced cross-section 60 mm in 12 measuring points *



* presented data relate to the measurement results of 12-02-2016
for cabinet BioTectum 1.2 with serial number 91512, Class II,
made by an accredited laboratory

The structure of the cabinet was created with highest quality materials with strict requirements for ergonomics, design and installation

- dtwo absolute HEPA filters with effectiveness 99.995 % for molecules $\geq 0.3\mu\text{m}$
- interior of the working area made entirely of stainless steel of class 0H18 (DIN 1.4301) in the seamless construction with rounded sides, which reduces to a minimum the number of surfaces presenting a risk of contamination
- V-shaped air inlet on the air curtain, which reduces the risk of blocking the curtain
- electrically actuated front glass, disposed at an angle relative to the worktop
- glassed sides of the cabinet
- armrest fixed along the whole length of working space, made of stainless steel
- bowl made of stainless steel put under the worktop, enables the collection of split liquids
- housing - electrostatically painted steel sheet in RAL9003 colour
- stylistic elements - orange colour RAL2000
- connections: combustible gas, inert gas or vacuum, power supply are placed on the top of the cabinet so there are no wires protruding beyond the contour of the cabinet. This design allows maximum use of laboratory space
- height of the working hole of the cabinet - 200 mm



Technical data

sound pressure level PN-EN ISO 11201	47,2 dB**
number of motors, type	3 electronically commutated motors
illuminance	800-1700 lx
mechanical vibrations PN-EN ISO 5349	≤ 0,005 mm
interior of the working area	stainless steel 0H18 (DIN 1.4301)
main filter and exhaust filter	type: H14 E ≥ 99,995 % according to EN 1822:2009 norm

Electrical specification

power supply	1N 230 V / 50 Hz
electrical sockets protection	10 A
power consumption in safe work mode	119 W
power consumption in SUSPEND mode	68 W

External dimensions

width	1286 mm
depth	795 mm
height with the frame	2100 mm + 30 mm


Dimensions of the working area

width	1200 mm
depth	600 mm
height	710 mm

** the measurement was taken using the following parameters: laminar air flow velocity 0,25 m/s; air flow velocity in the barrier 0,4 m/s; height of the cabinet hole 160 mm



*For your safety
the construction
of the cabinet was
designed to easily
clean both sides of
the glass closing the
working area*



- possibility of automatic and semi-automatic UV sterilization
- possibility of setting the sterilization with time delay from 0... to 1440 minutes
- UV sterilization programming in a weekly cycle
- programmable sterilization time: from 1...to 360 minutes



AutoProtect mode

Connecting the cabinet with emergency power generator UPS and activation of AutoProtect program allows you to continue work with the cabinet after power supply failure and gives the necessary time to properly secure the material. The power supply to the electrical sockets and solenoid valves is cut off in a split second whilst the proper, safe working conditions are maintained.

The program works in two variants:

1. Alerts the operator in case of power failure (default settings) and leaves the decision to continue work with the cabinet to the operator, including the decision to close the glass and turn off the unit.

2. Automatic mode - works according to the predefined settings:

- working time on emergency power supply UPS
- time after which the glass is automatically closed and the unit is turned off
- each step is signaled by an appropriate alarm message





We share the chemistry...

NATIONWIDE LABORATORY SUPPLY NETWORK

Alchem Grupa Sp. z o.o

CENTRAL REGION
warszawa@alchem.com.pl
Warszawa, ul. Sokratesa 17A
tel. 22 632 99 54

SUBCARPATHIA REGION
rzeszow@alchem.com.pl
Rzeszów, ul. Boja Żeleńskiego 25B
tel. 17 850 50 60

PODLASIE-WARMIA REGION
bialystok@alchem.com.pl
Białystok, ul. Włociańska 58
tel. 85 653 87 65

KUYAVIAN - POMERANIAN REGION
torun@alchem.com.pl
Toruń, ul. Polna 21
tel. 56 658 90 40

LOWER SILESIAN REGION
wroclaw@alchem.com.pl
Wrocław, ul. Giżycka 62
tel. 71 372 66 39

WEST POMERANIAN REGION
stargard@alchem.com.pl
Stargard, ul. 1-ej Brygady 35
tel. 91 573 00 34

GREATER POLAND REGION
wielkopolska@alchem.com.pl
Ostrów Wielkopolski, ul. Długa 19
tel. 62 736 80 81

SILESIAN REGION
bielskobiala@alchem.com.pl
Bielsko-Biała, ul. Kasprowicz 9/2
tel. 33 828 55 00

COMPANY'S HEADQUARTERS
alchem@alchem.com.pl
Toruń, ul. Polna 21
tel. 56 655 05 49